

# Count sheet

## T2 Alpha Femur Retrograde Indication Tray

Insert: 2355-3920

Ref. #	Description	Location	Qty.	Check bo
2351-0070	Tissue protection sleeve, long	Tray	1	
2351-0330	Screwdriver bit for condyles	Tray	2	
1806-0257	Revision Screwdriver Bit, T2 Femur	Tray	1	
2351-0140	Quick lock delta handle, modified trinkle	Tray	1	
2351-6312	Opening reamer trocar 11.5mm	Tray	1	
2351-6212	Opening reamer sleeve 11.5mm	Tray	1	
2351-0335	Countersink bit for condyles*	Tray	1	
2351-6000	Opening reamer handle	Tray	1	
2351-6112	Opening reamer 11.5mm	Tray	1	
2351-0145	Quick lock T-handle, modified trinkle	Tray	1	

### Bottom layer

Ref. #	Description	Location	Qty.	Check bo
2353-3903	Nail adapter femur retrograde	One tray level	1	
2353-3904	Nail holding screw femur retrograde	One tray level	2	
2353-3907	DTD femur retrograde	One tray level	1	
2351-0250	External compression device femur	One tray level	1	
0152-0218S	K-wire 1.8 x 310mm	One tray level	2	
2353-3908	PTD short femur retrograde	One tray level	1	
2353-3909	Knob for PTD short	One tray level	1	
2355-5030	Drills and pins insert	One tray level	1	

\*Available sterile and non-sterile; add an S to the ref# for sterile.

\* Please note that the instrument bank and implant set are separate part numbers and must be ordered separately.

Hospital signature \_\_\_\_\_

